Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25
32657 Lemgo - Germany
Telefon +49 (0)5261 701-0
Telefax +49 (0)5261 701-289
info@kometdental.de
www.kometdental.de

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EN

Mechanical tapered instruments for endodontic treatments (root canal files), type 2, in compliance with ISO 3630

1. Intended use

The instruments are used during endodontic treatments for the mechanical removal of diseased tissue and root filling material. Root canal treatments are indicated for example in case of an acute pulpitis due to profound caries and periapical lesions. Contra-indications include the preparation of milk teeth, abrupt apical curvatures, s-shaped curvatures and pronounced calcification of the root canal. The Openers / FQ08 are not authorized for use in the apical third of the root canal. In addition, the RE10 must not be used in the central third of the root canal, either, in order to prevent iatrogenic perforations. The instruments are used for the treatment and alleviation of the above listed diseases in combination with a dental power system in a rotary or reciprocating mode. They are intended for use by medical professionals for removing substance in patients with permanent teeth in a dental practice/office or clinic.

2. Warnings

- Irrespective of the file system used, the following side effects and incidents can occur at any point during endodontic treatments: Fracture of the instrument inside the root canal, vertical fracture of the root, periodontal endodontic lesions, open apex or perforation.
- Possible residual risk of injury to the patient and the user caused by the sharp and pointed instruments. The user has to take appropriate protective measures.
- Due to the delicate shape of the instruments, there is a risk of file fracture. This can be minimized by following the below listed instructions:
- Observe the recommendations concerning the use of mechanical NITI instruments and the cleaning and disinfection of the root canal.
- The cleaning and disinfection of the instruments has to take place mechanically in a cleaner/disinfector. Manual methods, including the use of an ultrasonic bath, should only be applied if this is authorized by the legal prescriptions valid in your country and if no automated procedure is available. Manual reprocessing is considerably less effective and reproducible than mechanical methods.
- Multiple use instruments may only be reprocessed 10 times at most.
- The maximum number of times a file can be used / number of root canals that can be prepared with one file indicated in the above table must not be exceeded. Ignoring this recommendation can reduce the efficiency or lead to fracture of the instruments.
- Clean the instruments mechanically during the treatment, for example in the interim bur block 595, and rinse the root canal thoroughly in order to remove any debris.
- For root canals with pronounced curves (>30°), narrow radius or multiple curves, the maximum permissible treatment is 1 root canal.

- Check the instruments during the treatment, i.e. after intermediate cleaning, for plastic deformation, such as untwisted spirals on the working part. Discard damaged instruments if applicable.
- It is recommended to define the canal morphology by means of x-ray imaging.
- Put the instruments in operation before inserting them into the canal. The instruments are used in a "pecking motion", i.e. in short up and down movements.
- The instruments may only be used by professional health care staff on human patients with permanent dentition at a dental practice or clinic.

3. Appropriate use

To ensure a successful treatment, the predefined treatment flow has to be followed and the specifications listed in the above table, such as the max. speed and the max. number of uses etc., have to be observed. Rotary root canal files are intended for use in a contra-angle with torque limit (see table above). Reciprocating root canal files are intended for use in established left-cutting reciprocating power units (see table above). Thanks to the special properties of the FQ (FQ08, FQ06, FQ04) and Procodile Q files, the instruments are highly flexible to facilitate pre-bending of the working part. As a result, the instruments might already be slightly bent when delivered. This does not impair the quality or proper function of the files

1. Exposure of the canal entrance with the Opener
2. Clearly expose the canal entrance, for example with a suitable opener (FQ08, OP, OPR, AK10), intended for use only in the coronal third of the root canal. In case of gently curved canals, the openers OP / FQ08 and OPR can be used up to the center third of the root canal.
3. Probing with manual files
4. Probing of the canal with manual root canal files, such as K-Files.
5. Mechanical creation of a glide path with the PathGlider
6. Optional: Mechanical creation of a glide path, e.g. with the PathGlider: PG03 / AK03 / FQ03 in size 015 requires previous probing with a manual file in size 010
7. PG03 in size 020 requires previous probing with a manual file in size 015
8. Mechanical preparation of the root canal
9. Instruments for the preparation of the root canal (FQ, F360, F6 SkyTaper, Procodile Q, Procodile) are used according to the single-length principle (long the entire working length).
10. Generally speaking, the size of the first Komet preparation file depends on the size of the probing file. The rule is to skip one ISO size. Example:
 - Probing /glide path with ISO 010: Use of a preparation file in size 020
 - Probing /glide path with ISO 015: Use of a preparation file in size 025 etc.
11. Retreatment:
 - Use of a solvent, if required, for example eucalyptus oil to soften the gutta-percha
 - Endo ReStart Opener, to be used in the coronal third of the root canal only, to remove the root filling. The instrument has a cutting tip with which it can even penetrate hard root canal fillings. The toothbrush is short because the use of the instrument is limited to the coronal third of the root canal.
 - Endo ReStart files, removal of the residual root filling according to the single-length principle. The instruments are provided with non-

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cutting tips that have a scraping function instead which enables them to penetrate even hard root fillings. The toothbrush is provided with pronounced turns to enable the file to screw itself into the root canal. To guarantee a sufficiently long service life of the Endo ReStart files, the gutta-percha has to be removed gradually and the file may only be used applying very low axial forces.

4. Reprocessing

- > Sterile packed instruments for single use
- The instruments are supplied in sterile packaging and can be used without prior processing, provided that the packaging is intact and without damage. Do not use the instruments if the packaging is damaged. Single-use articles (marked with a crossed-out 2 on the packaging) are not authorized for reuse. A risk-free reuse of the instruments cannot be guaranteed due to the risk of infection and/ or because the safety of the products can no longer be assumed.
- > Unsterile packed instruments for single use
- Those instruments that are supplied non-sterile have to be prepared prior to first use. Follow the procedure described under „reprocessing steps“. Silicone stoppers are intended for single use and have to be removed from the files before these are reprocessed. Single-use articles (marked with a crossed-out 2 on the packaging) are not authorized for reprocessing. A risk-free reuse of the instruments cannot be guaranteed due to the risk of infection and/ or because the safety of the products can no longer be assumed.
- > Sterile packed instruments for repeated use
- The instruments are supplied in sterile packaging and can be used without prior processing, provided that the packaging is intact and without damage. Do not use the instruments if the packaging is damaged. Follow the procedure described under „reprocessing steps“. Silicone stoppers are intended for single use and have to be removed from the files before these are reprocessed and be replaced by commercially available stoppers before the next use.
- > Unsterile packed instruments for repeated use
- Those instruments that are supplied non-sterile have to be prepared prior to the first and after every subsequent use. Follow the procedure described under „reprocessing steps“ below:
- Reprocessing steps:
 - Initial processing at the site of use and preparation: After use, place the instrument in a disinfectant/disinfection tank filled with a disinfecting and cleaning solution (e.g. DC1, co. Komet [validated procedure: 3 % solution, 5 minutes]) or Cidex OPA, co. Johnson & Johnson [validated procedure: undiluted]. Consult the instructions issued by the manufacturers of these agents for directions of use (concentration, temperature and max. and min. immersion time). Reprocessing should take place shortly after use (≤ 2 hours). The instruments have to be in the cleaning/disinfection tank during transport to the reprocessing site. Remove the instruments from the cleaning/disinfection tank immediately before starting the cleaning program. Remove any residual contamination and any residues from the cleaning/disinfection agent with the help of a nylon brush by rinsing the instrument for 1 minute under running water (less than 45 °C to prevent protein fixation).
 - Cleaning /disinfection: Place the instrument inside a suitable bur block in the washer /disinfector (EN ISO 15883) approved by DGHM or FDA) in a horizontal position, so that the spray jet directly hits the instruments through the openings. Start the Vario TD program, washer /disinfector, co. Miele or universal program, co. Melag. Thermal disinfection takes place in a washer /disinfector (at least 5 minutes at 90 °C or A₁ value ≥ 3000). Observe the instructions provided by the manufacturer with regard to the concentration and quantity of the detergent (e.g. Neodisher MediClean forte, co. Dr. Weigert). After the complete cycle, allow the instruments to cool and remove from the washer /disinfector.
 - Drying /control: Dry instruments with low-germ, oil-free compressed air, insisting on hard-to-reach areas. The visual inspection is done with the naked eye without magnification. In case of visible residual contamination, repeat the steps pre-treatment, cleaning/disinfection and drying. Do not reuse any instruments that show signs of damage, chipping, deformation or corrosion.
 - Sterilization: Pack the cleaned, disinfected and dried instruments into suitable disposable sterilization packaging. The sterilization has to be carried out in the autoclave type B (EN 13060/ EN 285 or ANSI AAMI S179/ approved by the FDA) at 134 °C with a hold time of at least 5 minutes and a subsequent drying time: 20 minutes. Devices that function according to a different process are less efficient and should therefore not be used. If the user applies a different method, the user has to validate the efficiency of such a method.

Please also observe the legal provisions valid in your country.

The efficiency of the described procedure was confirmed by an accredited laboratory. In case of any other procedure applied by the user, the efficiency needs to be validated by the user. The operator of medical products is responsible for ensuring that the reprocessing procedure is carried out with suitable equipment and suitable material by qualified personnel and in compliance with any legal provision valid in your country. Information on proper reprocessing of critical instruments according to EN ISO 17664, the present instructions for use, as amended, and more detailed information can be downloaded from our website <https://fr.kometdental.com> or requested from the manufacturer Gebr. Brasseler.

5. Storage

Store wrapped devices protected from UV radiation and high temperatures in a clean and dry environment. Do not store in the same room as solvents or chemicals.

6. Disposal

The devices have to be disposed of in tight, non-breakable and puncture-proof containers (protection from contamination).

7. Safety and liability

Please also observe the General Instructions for Use and Safety Recommendations in the current Komet Dental catalogue. The user is obliged to check the devices on his own responsibility prior to use to ensure that they are suited for the intended purpose. The correct application of the devices is the responsibility of the operator and user. In case of contributory negligence by the user, Gebr. Brasseler GmbH & Co. KG partially or totally declines liability towards the user for all resulting damages, particularly if these are due to the non-observance of the recommendations for use or other misuse by the user.

Store out of children's reach. For dental use only.

8. Procedure in case of serious adverse events

In case of serious adverse events, please proceed in line with the legal regulations in your country.

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© Gebr. Brasseler GmbH & Co. KG